
Full DNA Adapter Kit for MGI

RK21676 MiniSet (8 indices)

RK21677 MidiSet (24 indices)

RK21678 Set_A (48 indices)

RK21679 Set_B (48 indices)

www.abclonal.com.cn

version: N16I17v4.0

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Full DNA Adapters Kit for MGI MiniSet (RK21676)

组分名称	Barcode 序列	组分名称	Barcode 序列
MGI Full Adapter 1	TAGGTCCGAT	MGI Full Adapter 13	CGGCAATCCG
MGI Full Adapter 2	GGACGGAATC	MGI Full Adapter 14	ATCAGGATTC
MGI Full Adapter 3	CTTACTGCCG	MGI Full Adapter 15	TCATTCCAGA
MGI Full Adapter 4	ACCTAATTGA	MGI Full Adapter 16	GATGCTGGAT

Full DNA Adapters Kit for MGI MidiSet (RK21677)

组分名称	Barcode 序列	组分名称	Barcode 序列
MGI Full Adapter 1	TAGGTCCGAT	MGI Full Adapter 45	TGGACTTGGC
MGI Full Adapter 2	GGACGGAATC	MGI Full Adapter 46	GCTTGACAGG
MGI Full Adapter 3	CTTACTGCCG	MGI Full Adapter 47	AAGACCTCTA
MGI Full Adapter 4	ACCTAATTGA	MGI Full Adapter 48	AGTTGCCATA
MGI Full Adapter 13	CGGCAATCCG	MGI Full Adapter 57	ATCAACGGA
MGI Full Adapter 14	ATCAGGATTC	MGI Full Adapter 58	AACTGTACTG
MGI Full Adapter 15	TCATTCCAGA	MGI Full Adapter 59	GTACCTCAAT
MGI Full Adapter 16	GATGCTGGAT	MGI Full Adapter 60	GACTTCTAAT
MGI Full Adapter 41	TTAGATGCAT	MGI Full Adapter 61	TGAAGCGTTG
MGI Full Adapter 42	GTCCAGAGCT	MGI Full Adapter 62	CGTGCGATCC
MGI Full Adapter 43	CACGTGATAG	MGI Full Adapter 63	TCGGAAGGCA
MGI Full Adapter 44	CCACTAGTCC	MGI Full Adapter 64	CCGATGTCGC

Full DNA Adapters Kit for MGI Set A (48 indices) (RK21678)

组分名称	Barcode 序列	组分名称	Barcode 序列
MGI Full Adapter 1	TAGGTCCGAT	MGI Full Adapter 65	ACTTAGAATG
MGI Full Adapter 2	GGACGGAATC	MGI Full Adapter 66	TCCAAGCCTG
MGI Full Adapter 3	CTTACTGCCG	MGI Full Adapter 67	AGACGATGAT
MGI Full Adapter 4	ACCTAATTGA	MGI Full Adapter 68	CTCACAAGAC
MGI Full Adapter 13	CGGCAATCCG	MGI Full Adapter 69	CGTTCCTACT
MGI Full Adapter 14	ATCAGGATTC	MGI Full Adapter 70	GTGGTTGTGA
MGI Full Adapter 15	TCATTCCAGA	MGI Full Adapter 71	GAAGGCCTGC
MGI Full Adapter 16	GATGCTGGAT	MGI Full Adapter 72	TAGCTTGCCA
MGI Full Adapter 41	TTAGATGCAT	MGI Full Adapter 73	GACAATGCTC
MGI Full Adapter 42	GTCCAGAGCT	MGI Full Adapter 74	GCTAATCACA
MGI Full Adapter 43	CACGTGATAG	MGI Full Adapter 75	AGTCCATAGG
MGI Full Adapter 44	CCACTAGTCC	MGI Full Adapter 76	CTATCGCCTA
MGI Full Adapter 45	TGGACTTGGC	MGI Full Adapter 77	ATCGTGGTCT
MGI Full Adapter 46	GCTTGACAGG	MGI Full Adapter 78	TGGCTAATAC
MGI Full Adapter 47	AAGACCTCTA	MGI Full Adapter 79	CAGTGCAGAG
MGI Full Adapter 48	AGTTGCCATA	MGI Full Adapter 80	TCAGGCTGGT
MGI Full Adapter 57	ATTCAACGGA	MGI Full Adapter 81	ATACTCACGC
MGI Full Adapter 58	AACTGTACTG	MGI Full Adapter 82	ATGCTCCGCG
MGI Full Adapter 59	GTACCTCAAT	MGI Full Adapter 83	TGTGAAGTGTG
MGI Full Adapter 60	GACTTCTAAT	MGI Full Adapter 84	GAGAGGTGCT
MGI Full Adapter 61	TGAAGCGTTG	MGI Full Adapter 85	TGCACTGTAA
MGI Full Adapter 62	CGTGCATCC	MGI Full Adapter 86	GCCTAGGCAA
MGI Full Adapter 63	TCGGAAGGCA	MGI Full Adapter 87	CCATCATAGC
MGI Full Adapter 64	CCGATGTCGC	MGI Full Adapter 88	CATGGTAATT

Full DNA Adapters Kit for MGI Set B (48 indices) (RK21679)

组分名称	Barcode 序列	组分名称	Barcode 序列
MGI Full Adapter 89	CACCATGTCT	MGI Full Adapter 25	TAGAGGACAA
MGI Full Adapter 90	ATATGTCTGG	MGI Full Adapter 26	CCTAGCGAAT
MGI Full Adapter 91	AAGGAAGCGT	MGI Full Adapter 117	ATGTCTATCC
MGI Full Adapter 92	TCAAGACGTC	MGI Full Adapter 28	GCTGAGCTGT
MGI Full Adapter 93	CCGCTCAGTA	MGI Full Adapter 29	AACCTAGATA
MGI Full Adapter 94	GGTGTGTACA	MGI Full Adapter 30	TTGCCATCTC
MGI Full Adapter 95	TTCACGTAAG	MGI Full Adapter 114	CCAGAGTCAG
MGI Full Adapter 96	GGTTCCACAC	MGI Full Adapter 32	CGCTATCGGC
MGI Full Adapter 97	AGGTATTCTT	MGI Full Adapter 33	GCAACGATGG
MGI Full Adapter 98	CGAATGCAAC	MGI Full Adapter 34	TAATCGTTCA
MGI Full Adapter 99	TTCAACGGCG	MGI Full Adapter 35	GTTGCTCTA
MGI Full Adapter 100	CTCGGCGGAA	MGI Full Adapter 36	TCTCACACAT
MGI Full Adapter 101	ACGGTAATGG	MGI Full Adapter 37	CTGTTAGGAT
MGI Full Adapter 102	GATCCGACGT	MGI Full Adapter 38	CGCAGACGCG
MGI Full Adapter 103	TCACGATACA	MGI Full Adapter 39	AAGGATCATC
MGI Full Adapter 104	GATTCTCTTC	MGI Full Adapter 115	AACAGGCAGT
MGI Full Adapter 121	CCTTGATCAA	MGI Full Adapter 49	ATGTACGCAG
MGI Full Adapter 122	GGAAGTGGCA	MGI Full Adapter 50	TTAATGAGAT
MGI Full Adapter 123	AACATTCTAC	MGI Full Adapter 51	TGCGCCACTT
MGI Full Adapter 124	GACGCGAGTC	MGI Full Adapter 52	CATTAAGGCC
MGI Full Adapter 125	CTATAAACT	MGI Full Adapter 53	CCGCCTCAGA
MGI Full Adapter 126	AGTCTCGTGT	MGI Full Adapter 116	GCTCCATGAC
MGI Full Adapter 127	TCGGCCTATG	MGI Full Adapter 55	GCCGTTATC
MGI Full Adapter 128	TTGCAGACGG	MGI Full Adapter 56	GGAATATTGA

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